# FLIR camera test winter 2006

## Introduction

1. The objective is to record solar and lamp spectra at the same position as the Rockwell camera in a SPINOR setup to check for ability to flat field and sensitivity relative to the Rockwell.

### Setup

- 1. Connect the camera to a computer of choice, either Harry or a NSO machine. The camera connects to the computer via an Ethernet cable and a USB cable. Recommended Ethernet speed is 1000BT (Harry has 100BT).
- 2. Learn how to drive the FLIR application so that single 50msec exposures can be recorded onto disk in raw format.
- 3. Use the SPINOR setup for the Pietarila/Judge observing run
- 4. There is no need for the ASP computers, cameras, or modulator.
- 5. The SPINOR beam splitter needs to be installed behind the slit if not already
- 6. Adjust the 'L' fold mirror to put light on the grating if not already
- 7. Continue to use the 308.571/mm grating at the same incidence angle as used for the observing run, 56.35°.
- Remove the Rockwell camera from the HSG with great care. Remove and set aside the 1.5648μm pre-filter. Cover the camera entrance window with the translucent plastic hemispheric cover. Refer to the observing plan on the web for all the Rockwell items to pack.
- 9. Place the FLIR camera on a mount sideways so the long axis is vertical
- 10. Put the 1.56µm filter immediately in front of the entrance window of the camera. If it arrived with a lens, remove it.
- 11. Check that the slit width is  $50\mu m$
- 12. Adjust the camera position in focus until lines are sharp. It is OK to adjust the grating angle to center the 1.5648µm line.

#### **Measurements**

- 1. Record a couple of images on the Sun at 50msec exposure time.
- 2. Record a couple of dark images at the same exposure time.
- 3. Insert the calibration lamp, set the slit width to 1500 m and record a couple of images
- 4. Move the data to some place accessible

#### Comments

1. David Elmore 303 497 1580 (303 669 7264 mobile)